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Compliant with JP18 and USP <1207>

Packaging Integrity Testing of Pharmaceutical Products (CCIT)

Quantitative test methods that do not require wetting the product with water

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- Air leak testing device for sample inspections
- Barrier property checking (Oxygen, Steam, etc.)
- Establishing "Maximum allowable leakage limits"
- Testing is also quantified for scientific quality management
- Test Results that are both quantitative and reproducible
- Equivalent to a minimum opening diameter of $\phi 3 \mu\text{m}$
※Depending on measurement conditions
- GMP, DI, and Validation Compatible

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Leakage Testing Compliant with The Japanese Pharmacopoeia and USP <1207>

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Sealed Chamber Method
(Pressure Change Leak Test Method 1)

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Due to the improvement in the quality of pharmaceuticals and the internationalization of distribution, there is a growing demand for quantitative testing of packaging containers. JP18 includes references to "Integrity Assessment of Aseptic Pharmaceutical Packaging (G7-4-180)" and "Leak Test Methods for Aseptic Pharmaceutical Packaging (G7-5-180)", requiring quantification based on "maximum allowable leak limits," similar to USP <1207>. Both JP18 and USP <1207> require the use of positive and negative controls for integrity assessment. Fukuda recommends the use of pseudo-defects (such as pinholes and capillaries) with a specific pore size for the positive control.

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
When determining a "Maximum allowable leakage limit," use the MSP series and Artificial defect FPH series. Use cases include the following.

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Web site / Information page
There are download materials at the bottom.

https://fukuda-jp.com/column/packaging_container.html

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*1 "Test Parts" refers to the packaging containers being inspected, as well as samples. In this context, it is referred to as "Test Parts".

*2 If the test parts you want to test is larger than the size mentioned above, please consult us.

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Examples of Packaging Containers Compatible with Different Models and Detectable Pore Diameters

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Pseudo Defect FPH Series (Sold Separately)

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

This is a pseudo-defect (reference test piece) with specified hole diameter and hole length.

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Pinhole utilization case (can be used for creating difficult fine defects (positive control))

包装容器の密封試験を自動化、数値化、乾式化
The Transition to Automated, Quantified, and Dry Packaging Seal Tests

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Packaging containers are designed to protect their contents by preventing moisture and gases moving in and out, and microbial ingress. If the seal integrity is compromised with pinholes and any other defective sealing, quality will suffer due to degradation or degeneration of the contents. The MSP series provides a reproducible and quantifiable method, and conversion from the traditional water immersion tests.

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For air leak testing of packaging containers, leakage tests of the object to be tested are performed using pressure variations from the test parts. Large and small leaks are detected by measurement in the following two processes.

Testing Process for Large Leaks : Checking for large leaks such as coming away, peeling, and tearing of the packaging.

Testing Process for Small Leaks : Checking for defects (such as pinholes in film materials) that are often difficult to find by visual inspection.

Prior testing of the shape, volume, and test pressure, etc. of your test parts is required in order to perform optimal leak testing. Please do not hesitate to contact us to discuss your specific needs.

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Inspection Process One device can support two kinds of test methods.

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● **Data Integrity, Data Falsification Preventing (USB, Printer)**
Log data such as test results, pressure values, and setting values is encrypted, and can be recorded or converted to PDF format on a PC.

● **ID Control**
It is possible to identify the operator and the administrator separately using ID's and passwords.

● **Type Registration / Product Setting Function**
Product type and setting conditions can be registered in 1-100 (100 types) groups.

● **Self Check Function**
This function checks for malfunctions in various valves or sensor faults.

● **Barcode Function**
You can select a type, and input an ID at login.

● **Simple Setup**
Automatic settings of testing conditions using good / defective products.

● **Sensitivity Check**
Sensor sensitivity can be confirmed by built-in flow reference units (FFM series).

● **On-screen Display of Measurement Status**
Confirmation of measurement status according to data waveform, verification of setting value at time of start-up, and determination of the cause in the event of a problem occurring are available. Settings and operation are easy thanks to the touch panel screens.

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Number		Item	Description
①	0101	Model	Please check the specifications below for the test parts size.
	0102		
—	W	Measurement Method	Tank Pressurization Type/ Tank Partial Pressure Type
—	W	Pressure Range	Compound Pressure (Equipped with APU)
—	N	Additional Tank	None
②	V	Validation	With Validation
	N		Without Validation

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

Item	Description		
Measurement Method	Tank Pressurization Type/ Tank Partial Pressure Type Sealed Chamber Method (Pressure Change Leak Test Method 1)		
Test Pressure	1 ～ 200 kPa abs Compound Pressure (Equipped with APU)		
Products Setting	1 ～ 100 Group (100 Types)		
Operation Temperature/ Humidity Range	0 ～ 40 ℃ 45 ～ 85 % RH (Without condensation)		
Power Source/ Applied Current	100 ～ 240 VAC 50/ 60 Hz 100 VA		
External Input/ Output	Serial Communication 3 Port (PC, Printer, Barcode reader) USB 1 Port Ethernet 1 Port		
External Dimensions/ Mass	MSP-0101	W380 × D560 × H290 mm	25 kg
	MSP-0102	W380 × D560 × H470 mm	30 kg
	Dedicated Cradle (Sold Separately)	W380 × D560 × H700 mm	55 kg
Test Parts Size	MSP-0101	W130 × D180 × H65 mm or below	
	MSP-0102	D φ 120 × H220 mm or below	
※ Please consult FUKUDA if testing is required for sizes larger than that shown here.			
Pneumatic Pressure Source	Pneumatic Pressure Source 0.4 MPa or above. Please supply clean, stable, and dry air. Recommended conditions according to ISO 8573-1: 2010. Compressed air purity classes 1, 3, 1		
Utility Dedicated Cradle (Sold Separately)	Air leak testing is test method using compressed air. If the condition of dry air used is low, the test result will be affected. In order to ensure testing accuracy, we provide a dedicated cradle that meets the recommended pneumatic source conditions. The dedicated cradle is made up of a compressor, vacuum pump, and a dryer.		
Optional (Sold Separately)	● FUKUDA Pinhole – FPH series (φ 5,10,15,20 μm etc.) ● Pneumatic unit FRZ-0057 with Bourdon tube ● Vacuum pump ● Barcode reader ● USB memory ● Printer (thermal/ dot impact type) ● Printer roll paper		

● Options such as a core, experiment fee, utility dedicated cradle (vacuum pump, compressor, etc.) are not included in the price of the unit.

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● **Total inspection system**
For inquiries regarding the full inspection equipment system, please contact our website or the specialized sales team for pharmaceutical packaging.

●写真は参考例です。The photo shows an example for reference purposes.